# OC SERIES OXYGEN CONCENTRATOR

Edition: V1.0

User Manual

SYSMED(CHINA) CO.,LTD  
11-2-3, No.17 WenSu Street, HunNan New District,  
ShenYang, China

**GCE**<sup>®</sup>  
Healthcare

# OC SERIES

English .....	4	EN
Deutsch .....	26	DE
Français .....	50	FR
Nederlands .....	74	NL
Español .....	98	ES
Português .....	122	PT
Italiano .....	146	IT
Čeština .....	170	CS

# ENGLISH

## User manual: M50

SPECIAL NOTES .....	5
WARNINGS.....	5
CAUTIONS.....	7
RADIO FREQUENCY INTERFERENCE.....	7
DANGER.....	7
ACCESSORIES.....	7
1. Product overview .....	8
Working principle .....	8
Flow diagram.....	8
Application scope.....	8
2. Unpackig, transport and storage .....	8
Unpacking .....	8
Transport and storage .....	9
Using conditions.....	9
3. Technical specification.....	9
Model and function definition.....	9
4. Part name & function .....	10
5. Operation instructions .....	13
Operation instructions .....	13
Button board definition .....	14
Display instruction .....	14
Additional optial function.....	17
6. Maintenance .....	18
Routine maintenance .....	18
Professional maintenance .....	19
7. EMC information.....	19
8. Trouble shooting guide .....	22
9. Quality warranty .....	23
Attachment I circuit diagram.....	24
10. Glossary .....	25

## SPECIAL NOTES

 **DEALER:** This manual **MUST** be given to the end user of this product.

 **USER:** Before using this product, read this manual carefully and save for future reference. Please read this manual carefully, **DO NOT** use this product or any available optional equipment without first completely reading and understanding these instructions and any additional instructional material supplied with this product or optional equipment, otherwise, injury or damage may occur.

If you are unable to understand the warnings, cautions or instructions, contact dealer before attempting to use this equipment.

- Consult with your physician before using the equipment. Use this concentrator as prescribed by your physician. Do not change the settings of flow and duration before consulting your physician to avoid insufficient oxygen supply or carbon dioxide retention.
- This equipment is not to be used as life supporting or life sustaining device. The device can only be use for supplemental oxygen. It's necessary to have a backup oxygen supply (such as oxygen cylinder or liquid oxygen) in the event of a mechanical failure or power outage.
- Close supervision is necessary when this product is used on children or physically challenged individuals. Additional monitoring or attention may be required for patients using this device who are unable to hear or see alarms or communicate discomfort.
- Contact your dealer or physician immediately if any untoward effect occurred.
- If the concentrator is not working properly or if water spills into the machine, switch off the device and call your dealer for examination/ repair. Do not dismantle the device.
- Patients with severe carbon monoxide poisoning are prohibited to use the equipment.

## WARNINGS

 **WARNING:** A warning represents the possibility of harm to the operator or patient.

 **WARNING:** For proper operation, your concentrator requires unobstructed ventilation. The ventilation ports are located at the rear base of the device and at the side air inlet filter. Keep the device at least 30 cm away from walls, furniture, and especially curtains that could impede adequate airflow to the device. Do not place the concentrator in a small closed space (such as a closet).

 **WARNING:** Do not remove the covers of this device. Servicing must be referred to an authorized and trained SYSMED home care provider.

 **WARNING:** In the event of an equipment alarm or if you are experiencing any signs of discomfort consult your home care provider and/or your health care professional immediately.

 **WARNING:** Oxygen generated by this concentrator is supplemental and should not be considered life supporting or life sustaining. In certain circumstances oxygen therapy can be hazardous; any user should seek medical advice prior to using this device.

 **WARNING:** Where the prescribing healthcare professional has determined that an interruption in the supply of oxygen, for any reason, may have serious consequences to the user, an alternate source of oxygen should be available for immediate use.

 **WARNING:** Oxygen vigorously accelerates combustion and should be kept away from heat or open flame. This device is not suitable for use in the presence of a flammable anesthetic mixture with air, oxygen or nitrous oxide.

 **WARNING:** Do not smoke, allow others to smoke or have open flames near the concentrator when it is in use.

 **WARNING:** Do not use oil or grease on the concentrator or its components as these substances, when combined with oxygen, can greatly increase the potential for a fire hazard and personal injury.

 **WARNING:** Do not use the oxygen concentrator if either the plug or power cord is damaged. Do not use extension cords or electrical adapters.

 **WARNING:** Do not attempt to clean the concentrator while it is plugged into an electrical outlet.concentrator.

-  **WARNING:** Device operation above or outside of the voltage, flow, temperature, humidity and/or altitude values specified may decrease oxygen concentration levels.
-  **WARNING:** Your home care provider is responsible for performing appropriate preventive maintenance at the intervals recommended by the device manufacturer.
-  **WARNING:** The specified operating environment temperature range of this device is 10-37°C. Below 10°C, the compressor may be difficult to start, and above 37°C, the compressor may overheat, thus reducing device life. It takes a minimum 2 hours for the device to warm from the minimum storage temperature until it is ready for intended use, and 2 hours minimum for the device to cool from the maximum storage temperature until it is ready for intended use.
-  **WARNING:** Long oxygen tubing should be kept away from infants and small children to avoid entanglement which can cause suffocation due to strangling.
-  **WARNING:** For any liquid spilled on the device, turn off power and unplug from electrical outlet before cleaning up spill. "I" means Power on, "O" means Power off.
-  **WARNING:** Keep the device upright during use or transportation. Keep the device away from high pollution or high smog areas.
-  **WARNING:** Oxygen therapy is harmful under certain conditions. Sysmed recommends that you consult a doctor before using the device. Oxygen therapy should be regularly evaluated and the oxygen therapy settings set according to doctor's recommendations. Do not adjust oxygen flow at will to avoid oxygen deficiency, oxygen poisoning or carbon dioxide retention.
-  **WARNING:** The oxygen concentrator should be kept away from heat or spark sources (including static electricity) and open flames.
-  **WARNING:** Smoking is strictly prohibited during oxygen inhalation. Do not smoke in the same room as the oxygen concentrator or any oxygen-carrying accessories. If you want to smoke, you must turn off the oxygen concentrator and remove the cannula or tubing . Leave the room where the oxygen tube or mask or oxygen concentrator is located. If you cannot leave the room, you must turn off the oxygen concentrator for 10 minutes before smoking.
-  **WARNING:** Keep power cord away from heat or spark sources (including static electricity) and open flames. When not in use turn off the power to the device.
-  **WARNING:** DO NOT put oxygen tubing under a bedspread or chair cushion.
-  **WARNING:** Before cleaning the oxygen housing or replacing the fuse, the power plug must be unplugged to avoid electric shock.
-  **WARNING:** Geriatrics, children, or other patients who cannot communicate fluently may require additional monitoring to avoid harm. Excessive kinks in the nasal cannula may affect the oxygen output.
-  **WARNING:** Reserve an oxygen bottle at home for use when the oxygen concentrator does not work during a power outage.
-  **WARNING:** Oxygen tubing should be equipped with a fire safe valve which can interrupt the flow of oxygen in the event of a fire.
-  **WARNING:** Each SYSMED oxygen concentrator should be used with crush-proof, certified oxygen cannula and tubing fitted with a fire safe valve (supplied by the home oxygen provider). Tube length should not exceed 4 meters, and the fire safe valve should not be taken off during usage.

## CAUTIONS

A caution represents the possibility of damage to the equipment.

- Do not place liquids on or near the device.
- If liquid is spilled on the device, turn the power off and unplug from electrical outlet before attempting to clean up spill. Call your home care provider if device does not continue to work properly.
- Set the oxygen flow rate according to your doctor's prescription and adjust the float to the middle of the flow meter to the scale value, which is the set flow rate.

## RADIO FREQUENCY INTERFERENCE

This equipment has been tested and found to comply with EMC limits specified by IEC/EN 60601-1-2. These limits are designed to provide a reasonable protection against electromagnetic interference in a typical medical installation.

EN

## DANGER

To Reduce the Risk of Burns, Electrocution, Fire or Injury to Persons, please pay attention to the followings:

DO NOT disassemble. Refer servicing to qualified service personnel. There are no user serviceable parts.

Avoid using while bathing. If continuous usage is required by the physician's prescription, the concentrator MUST be located in another room at least 3 meters from the bath.

DO NOT come in contact with the concentrator while wet. DO NOT place or store product where water or other liquid can spill into the concentrator .

DO NOT reach for product that has fallen into water. Unplug IMMEDIATELY and contact dealer.

Use this concentrator within its technical specification as stated in the user manual.

 **WARNING: DO NOT connect the concentrator in parallel or series with other oxygen concentrators or oxygen therapy devices.**

## ACCESSORIES

SysMed products are specifically designed and manufactured for use in conjunction with SysMed accessories. Accessories designed by other manufacturers have not been tested by SysMed and are not recommended for use with SysMed products.

It may reduce the performance of the concentrator when use unspecified humidifier or other accessories.

## 1. PRODUCT OVERVIEW


### WORKING PRINCIPLE

Oxygen Concentrator for Medical Use adopts the principle of pressure swing adsorption technology. At normal temperature, the machine can continuously separate medical oxygen with high concentration (93%±3% oxygen for medical use) from the air when its power is turned on. The machine is easy to operate and quick to use, its flow can be adjusted. The oxygen supplying method is unique, and is superior to liquefied and high-pressure oxygen.

EN

**Note:** Using oxygen concentrator doesn't influence oxygen content in the air of the room.

### FLOW DIAGRAM

### APPLICATION SCOPE

This equipment is used to provide oxygen only for medical use and health care.

It is not be used for life supporting or life sustaining. SysMed recommends an alternate oxygen supply device in the event of power outage, alarm condition, mechanical failure, or for the persons who are in dire need of oxygen. (Such as oxygen cylinder or liquid oxygen).

**Contraindication:** Patients with severe carbon monoxide poisoning are prohibited to use the equipment.


### Packing Method

## 2. UNPACKING, TRANSPORT AND STORAGE

### UNPACKING

1. Check for any obvious damage to the carton or its contents. If damage is evident, notify the carrier, or your local dealer.
2. The carton is design into two parts (carton body and carton base), it can be easily unpacking by lifting the carton body upwards. Remove all foam filled block and loose packing from the carton.
3. Carefully take out all the components and the machine from the carton. Inspect/examine exterior of the oxygen concentrator for nicks, dents, scratches or other damages. Check all components as listed. If there are missing components or quality issue, please contact your local dealer.

**NOTE:** Please keep the carton and packing materials for storage and transport.

## TRANSPORT AND STORAGE

1. Keep the machine upright placement during transport, DO NOT place horizontally or upside down.
2. Ambient Temperature: -30°C-70°C.
3. Comparative Humidity: 15-95%RH(No condensation)
4. Atmospheric Pressure: 500 ~ 1060hpa

## USING CONDITIONS


- Electrical Requirements: 220V, 60Hz (115V, 50Hz – North America versions)
- Operating Environment Temperature: 10 ~ 37°C
- Relative Humidity: 20-65%RH
- Atmospheric Pressure: 700hPa ~1060hPa
- Environment: dry, well-ventilated, without dust, corruption and toxicological gas, away from sunlight and no intense electromagnetic interference.
- Keep unit at least 10cm away from walls, draperies, furniture, and the like.
- Altitude: Lower than 1828 meters, otherwise, it will reduce the oxygen purity.

**Note:** Add the voltage stabilizer if the power supply voltage fluctuates over -15%~+10%, before start the concentrator.

EN

## 3. TECHNICAL SPECIFICATION

### MODEL AND FUNCTION DEFINITION

### Form 1 function and parameters

PARAM- ETERS MODEL NO.	TEMPER		SOUND LEVEL (DB (A))	FUNCTION AND SETTING
	Flow (L/ min)	Purity (%)		
OC-E80	1-8	93±3%	≤50	Standard function: LCD display, oxygen concentration indicator (Red, green, yellow bar on the LCD), temperature alarm, high/ low pressure alarm, low oxygen concentration alarm, power failure alarm, timed setting, intelligent fault diagnosis, maintenance reminder. Additional optional function: <input type="checkbox"/> Positive pressure outlet
OC-E100	1-10	90±3%	≤50	


Parameter setting comparison form (For 220V machine)

MODEL NO.	OC-E80	OC-E100
Power Consumption (W)	480	480
Rating current (A)	2.25	2.25
Electrical Requirements	220V, 50Hz	220V, 50Hz
Flow rate (L/min)	1-8	1-10
Outlet pressure (Mpa)	0.05±0.005	0.05±0.005
Equipment Class and Type	Class IIa	Class IIa
Nebulized Particle (optional)	Median Diameter≤5μm, Percent: 90%	Median Diameter≤5μm, Percent: 90%
Net weight (kg)	27 kg	27 kg
Dimension (mm)	381x347x689	381x347x689
EMC Compliance to	IEC 60601-1-2	IEC 60601-1-2

EN

Notice: The information contained in this manual is subject to change without notice.

#### 4. PART NAME & FUNCTION

## FRONT VIEW

FIG. 1

EN

### COMPONENTS:

Main components of the concentrator include: compressor, filter, molecular sieve set, controlling panel, flow meter and humidifier. Additional optional parts: positive pressure outlet. Accessories are not shown.

#### 1. Humidifying Bottle

It is also called Oxygen humidifier, and is used to humidify the Oxygen in order that the dry Oxygen doesn't irritate the throat and nasal mucosa, and prevent phlegm from drying.

#### 2. Connecting Nut

#### 3. Oxygen Outlet

#### 4. Humidifier Connector

#### 5. Oxygen Flow-meter

The float in side is a sign to show the flow rate outlet from this concentrator.

#### 6. Knob of the Flow-meter

It is also called flow-regulating valve and used to adjust and control Oxygen flow rate. Please do not turn the button abruptly, or the valve spool may be broken. Turn on the valve anticlockwise and turn out it clockwise.

#### 7. LCD Display

It shows running time, elapsed time, oxygen purity bar (type OC-E), and real time oxygen purity (type OC-P);

Malfunction code of high/low pressure alarm, temperature alarm and other troubles alarm; And maintenance reminder icon etc. Guiding the users use the machine in a more scientific and safer way. (See Fig. 2)


FIG. 2


- 7.1 Oxygen purity percentage (type OC-P)
- 7.2 High purity indicator
- 7.4 Low purity indicator
- 7.5 Malfunction alarm
- 7.6 Maintenance reminder
- 7.7 Timed setting, Running time and Malfunction code display
- 7.8 Elapsed time

**EN**

- 8. Power Switch
- 9. Adjusting Panel and Buttons
- 10. Positive Pressure Outlet (See Fig. 1)
- 11. Model label
- 12. Over current protector
- 13. Power cord socket

- 14. Inlet filter access panel
- 15. Cabinet Air intake filters (two sides)

## 5. OPERATION INSTRUCTIONS

### OPERATION INSTRUCTIONS

1. The oxygen concentrator should be placed on flat ground with clear, ventilated environment to avoid smoke, corruption and toxicological gas, away from sunlight, heaters; Keep the unit all sides at least 30cm away from walls, drapery, furniture, and the like.

#### NOTE:

- DO NOT put sundries, water or oil container and other coverings on the top of this concentrator.
- DO NOT place anything at bottom of this concentrator. Ensure the bottom exhaust is not obstructed during operating, otherwise the machine will over-heat and cannot work normally.
- DO NOT place the unit in a narrow area or in a closet during

2. Install the inlet filter into the concentrator as steps below:

- a. Remove the filter access panel,
- b. Insert the inlet filter, turn clockwise to install,
- c. Reinstall the access panel.

**NOTE:** The paper filter needs to replace when it is dirty.

3. Press the button of humidifier connector to remove all set of the humidifier. Turn the humidifying bottle clockwise to remove cap from bottle, pour some distilled water to the scale mark, and replace the humidifier cap on the bottle, and then securely tighten.

#### NOTE:

- Use distilled water or cold boiled water in the humidifier, and refresh everyday.
- Water level is better to be in the middle of max and min line.
- To unload, turn the bottle clockwise.

4. Insert the humidifier connector to re-connect the humidifier to the machine. Turn the knob of the flow meter anticlockwise to ensure that it is unblocked.

5. Connect the power, insert the plug of power line into the power socket of the machine, and the other end of the plug connects with indoor power socket, turn on the power switch.

#### NOTE:

Be sure to use safe and approved sockets, the plug should be inserted into the socket tightly, if extend power strip is needed, please choose the product with electrical safety certification

6. The way to check whether the humidifier is properly tightened, use your finger to gently block the outlet of humidifier about 20 seconds, if the float in the flowmeter falls to the bottom, remove your finger and listen for a whistling sound from bottle. If so, it indicates that the humidifier is tight, if not, it indicates the humidifier leaks. Take off the humidifier, replace and tighten the cap again. If there is still no sound after adjustment, please contact your local dealer.

7. Regulate the output flow rate, adjust to float such that flow mark is in the middle of the float.

#### NOTE:

- Adjust the knob of the flowmeter, clockwise to decrease, anticlockwise to increase.
- For oxygen therapy, it is very important for the patient to strictly follow physician's advice, take oxygen treatment according to the flow rate and time set by physician. Do not adjust flow rate and time at will.

- Do NOT turn the ball in the flowmeter up to the red line, otherwise, the purity of oxygen will reduce.
8. Connect the intake end of oxygen tube to the outlet of humidifier, then set the nasal cannula over patient's ears, insert the nasal tube into patient's nostrils for starting the oxygen treatment. If longer tube needed, connect the extension tube with outlet of humidifier and oxygen tube.
  9. Turn off the power, when the oxygen therapy is finished.

**NOTE:**

EN

- Unplug the power plug, if continuous use is not needed.
10. OC series oxygen concentrator has the function of timing, please refer to the Timed setting operation introduction.
  11. OC series oxygen concentrator has the optional function of nebulizer, please refer to the introduction of nebulizer function.
  12. OC series oxygen concentrator has the function of oxygen concentration monitoring, when the machine starts working, three indicator bars of oxygen concentration on the LCD will light at the same time for 3 seconds and then will go out at the same time. After 6 to 7 minutes, the LCD indicates the normal concentration.

**BUTTON BOARD DEFINITION**

**1. Button control to increase timing**

Press the button one time for 10min up

**2. Button control to decrease timing**

Press the button one time for 10min down

(20 seconds after time setting, it will start to counting mode automatically.)

**3. Alarm Silence Button**

Press this button to stop the alarm sound.

**DISPLAY INTRODUCTION**

**1. Start within 5min display status**

The oxygen purity indicator, hour, minute and separator flicker at the same time.

The LCD displaying as below:

Model OC-S has no oxygen indicator

Model OC-E

## 2. Continuous operation mode


Under this mode, the LCD indicates elapsed time and current running time as display below:

EN

## 3. Timed Setting mode:

Set the time by pressing the increase or decrease timer buttons, each press will adjust the timer upwards or downwards by 10 minutes. When setting finished, wait for 20seconds it will start to count the time automatically. The concentrator will stop working, till the time runs out.

## 4. Oxygen Parameters Display

At rated flow, stated concentration should be achieved after initial warm-up period (less than 30 minutes). High O2 purity, Low O2 purity, display as below:


If oxygen purity is higher than 85%, the high O2 purity bar will be shown, and the concentration level is shown.

## 5. Malfunction state display

### 1) High Pressure Alarm


While the working pressure is higher than specification , the equipment will alarm. The malfunction code is H:07, and the concentrator will stop. LCD display as below:

EN


### 2) Low Pressure Alarm

While the working pressure is lower than specification , the equipment will alarm. The malfunction code is H:06, and the concentrator will stop. LCD display as below:


### 3) High Temperature Alarm

If the temperature of compressor is higher than specification, the concentrator will stop, and the malfunction code is H:05. LCD display as below:

### 4) The fault of OSD Sensor circuit


The malfunction code is E:03, the concentrator will continue work. LCD display as below:

### 5) Routine Maintenance Reminder

When the machine's working time reach to the maintenance time, this concentrator will remind users when routine maintenance is due as, display as below

EN

After serviced by professional, the reminder alarm will stop.

### 6) Power Failure Alarm

It will alarm if the power failure while the concentrator is working.

## ADDITIONAL OPTIONAL FUNCTION

### Positive pressure outlet

The Oxygen concentrator provides nebulizer point, with optional simple nebulizer, it can be used for aerosol therapy. The inner diameter of the connect tube for nebulizer and the joint is  $\Phi 5-5.5\text{mm}$ . Nebulizer specification: NebEasy 3558-C02.

Operating steps as below:

1) As Fig.3, connect the tube with tube joint, insert tube joint to the nebulization outlet connector of oxygen concentrator.

FIG. 3

Tube joint

2) As Fig.4 & Fig.5, turn on the power switch, close the flow meter and connect the other side of the connection tube with the joint of medicine filled nebulization nozzle, and then nebulization treatment can begin.

FIG. 5

FIG. 4

3) As Fig.6, when the therapy is finished, take off the tube. Press down the external fixation steel rim, and take out the tube joint.

FIG. 6


EN

**Note:** Nebulizer must be cleaned after each use. Clean nebulizer and connection tube a soft cloth dampened with detergent and clean water; for the nebulization nozzle and mask use clean water to clean first, then disinfect and sterilize by dipping into medical alcohol for 5min, and clean water, and finally put them in the packet after drying. (Detailed clean method please read the simple nebulizer user manual)

## 6. MAINTENANCE

### ROUTINE MAINTENANCE

1. While powered off clean the outside case with soft damp towel and little detergent, and then wipe it up with dry towel, once or twice per month.
2. Air intake filter is used to prevent dust, and need to be cleaned regularly. It is a critical maintenance step to clean the filter, at least once a week depending on environmental conditions. Operation steps as below.
  - 1) Take off the two filter nets on both sides of the case, take out the sponge mesh.
  - 2) Rinse the sponge mesh with clean water. Get ride of extra water and dry naturally.
  - 3) Set back after drying for future use.
  - 4) Re-install the sponge the and the filter nets.

### NOTE:

- If replacement is needed, please use accessories supplied by Sysmed. Accessories designed by other manufacturers have not been tested by Sysmed and are not recommended for use with Sysmed products.
  - Do not operate the unit without the air intake filter in place.
3. Over-load protector
- Over-load is caused by over heating, after checking and trouble shooting, press “reset” to re-start the concentrator.

4. Cleaning the Humidifier and cannula
- 1) Clean the cannula every day. First clean with detergent and then rinse with clean water completely, and dry naturally. Changing a new cannula once a month is recommended.
  - 2) The distilled water or cold boiled water in humidifier should be refreshed every day. The humidifier should be cleaned and disinfected every 3 weeks, clean the humidifier with detergent and then rinse with warm or hot water, dip it in household disinfectant for 20-30 minutes, clean again with warm or hot water, and then dry naturally. If not use, put it in a clean bag.
  - 3) If oxygen tube or cannula is needed, please contact your local home care provider or dealer for recommendations. They should also give you instruction on the proper usage, maintenance, and cleaning.

## PROFESSIONAL MAINTENANCE SERVICE

You should ask for professional maintenance from appointed provider or dealer every 4,000 hours.

**NOTE:** Only the professional or dealer trained and authorized by the manufacturer can do the repairs or adjustment.

## 7. EMC INFORMATION

This product is classified in Class I, type B. Not for life supporting.

Manufacturer statement: the product should be used in the following electromagnetic environment; the user should ensure that the product should be used in the following electromagnetic environment.

EN

EMISSION TEST	CONFORMITY	ELECTROMAGNETIC ENVIRONMENT
Radio-frequency emission CISPR 11	Group 1	The device uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
Radio-frequency emission CISPR 11	Class B	The device is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network.
Harmonic emission IEC 61000-3-2	Class A	
Voltage fluctuations/ Flicker emissions IEC 61000-3-3	Complies	
Emission of Radio Frequency Energy RTCA/DO-160G Section 21	Category M	This device is suitable for use onboard commercial airplanes inside passenger cabin.

**Note:** Oxygen concentrator cannot be used with other machines in a short distance. If it's has to be used in such an environment, it has to be tested.

Manufacturer statement: the product should be used in the following electromagnetic environment; the user should ensure that the product is be used in the following electromagnetic environment.

IMMUNITY TEST	IEC60601 LEVEL	LEVEL	ELECTROMAGNETIC ENVIRONMENT
Electrostatic Discharge (ESD) IEC 61000-4-2	±15kV air ±8kV contact	±15kV air ±8kV contact	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical Fast Transient/Burst IEC 61000-4-4	±2kV for power supply lines ±1kV for input-output lines	±2 kV for supply mains ±1 kV for input/output lines	Mains power quality should be that of a typical home or hospital environment.

Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical home or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% UT (>95% dip in UT) for 0.5 cycle at 45 degree increments 70% UT (30% dip in UT) for 0.5 seconds <5% UT (>95% dip in UT) for 5 seconds	<5% UT (>95% dip in UT) for 0.5 cycle at 45 degree increments 70% UT (30% dip in UT) for 0.5 seconds <5% UT (>95% dip in UT) for 5 seconds	Mains power quality should be that of a typical home or hospital environment. If the user of the device requires continued operation during power mains interruptions, it is recommended that the device be powered from an uninterruptible power supply or a battery.
Power frequency(50/60Hz) magnetic field IEC 61000-4-8	30A/m	30A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical hospital or home environment.

NOTE: UT is the a.c. mains voltage prior to application of the test level.

Conducted Radio-frequency IEC 61000-4-6 Radiated Radio-frequency IEC 61000-4-3	3V 150kHz-80MHz 6 Vrms Amateur Radio & ISM Bands between 150 kHz and 80 MHz 10V/m 80MHz-2.7GHz	3V 150kHz-80MHz 6 Vrms Amateur Radio & ISM Bands between 150 kHz and 80 MHz 10V/m 80MHz-2.7GHz	Portable and mobile RF communications equipment should be used no closer to any part of the device, including cables. Suggested the distance: d=1.2 150kHz-80MHz d=1.2 80MHz-800MHz d=2.3 800MHz-2.5GHz In which P is according to the regulations of transmitter manufacturers transmitter maximum output power, in watts (W) as the unit, isolation distance d is recommended, in meters (m).b Fixed field strength of rf transmitter based on electromagnetic survey to determine a, should conform to the level that are lower than in each frequency range. Interference may occur in the vicinity of equipment marked with the following symbol: 
---	--	--	--

Note 1: in the 80 MHz and 80 MHz frequency, using high frequency formula.

Note 2: these guidelines may not be suitable for all situations, the electromagnetic propagation by buildings, objects and the effect of absorption and reflection of the body.

A: stationary transmitter field, such as wireless cellular/cordless phones and ground mobile radio base station, amateur radio, AM or FM radio and television broadcasting, etc., the field intensity in theory can predict. For the evaluation of fixed RF transmitter electromagnetic environment, electromagnetic site survey should be considered. If the equipment can be measured field intensity is higher than the place of the application of radio frequency (RF) in line with the level, the machine should be observed to verify their works. If the observed abnormal performance, supplementary measures may be required, such as to directional or location on the machine.

B: in all 150 KHZ ~ 80 MHz frequency range, the field strength should be lower than 30 v/m

EN

Portable and mobile radio communication equipment and equipment or system isolation distance between recommendation: the equipment is expected in the RF radiation harassment is controlled under the electromagnetic environment of use. Based on maximum output power of communication equipment and the equipment of the buyers or users with the following recommended keeping a portable and mobile radio communication equipment and the minimum distance between the machines to prevent electromagnetic interference.

TRANSMITTER RATED POWER/W	DIFFERENT DISTANCE/M		
	150KHZ~80MHZ $D=1.2\sqrt{P}$	80MHZ-800MHZ $D=1.2\sqrt{P}$	800MHZ-2.5GHZ $D=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For not listed above rated maximum output power transmitter, recommended isolation distance d, in meters (m), can use the corresponding column in the formula to determine the transmitter frequency, P here is provided by the transmitter manufacturers transmitter, maximum output power rating in watts (W) as the unit.

Note 1: in 80 MHz and 800 MHz, formula is used at a high frequency range.

Note 2: these guidelines may not be suitable for all situations, the electromagnetic propagation by buildings, objects and the effect of absorption and reflection of the body.

#### POWER CORD FOR OXYGEN CONCENTRATOR

Whether shielding	No	
Type	125V, 10A	250V, 16A
Length	1.8 m	

Warning: Using the specified accessories and power cord may lead to the increase of oxygen concentrator launch or low immunity. Other cables and accessories may negatively affect EMC performance.

Basic performance:

1. Oxygen concentration: before/after test, the oxygen concentration is qualified.
2. Oxygen flow rate: before/after test, the oxygen flow rate is qualified.
3. Monitor function: work normal, qualified.

## 8. TROUBLE SHOOTING GUIDE

EN

SYMPTOM	PROBABLE CAUSE	SOLUTION
Power on, the equipment is not working but LCD display well	Start capacity of compressor is broken or compressor is not working	Call service provider or dealer
Power on, the equipment is not working, LCD does not display, or works intermittently	Power cord may not be plugged in properly or bad contact	* Check whether the power cord is damaged or not. * check whether the power cord is plug well. If no above question, call service provider
Exhaust sound too loud	* Outlet muffler taken off * Outlet muffler broken	Call service provider
No exhaust sound but safety valve has air out	* Pneumatic valve broken * controlling panel broken	Call service provider
No oxygen outlet or the outlet flow is too small	* oxygen cannula kinked or blocked * humidifier bottle not tightened	* unknit the cannula * re-install the humidifier cap Call service provider, if still cannot solve the trouble
The ball in flow meter is uncontrolled by the knob	* the flow knob is not tightened * turn the knob abruptly or too fast	* to tight the knob,turn the knob slowly and softly Call service provider, if still cannot solve the problem
Water back to cannula	* temperature difference caused by the weather, too near to the wall, drapey or furniture, and the location of the equipment and cannula is different	* dry the inside of humidifier cap * Do NOT use hot water * Do NOT over-filled humidifier * keep the same temperature of equipment and cannula (at same room)
Over-heating	* cooling fan is not working * compressor is broken	* put your hand on the side of inlet, to feel whether has cooling wind enter, yes means cooling fan is well, no means it's broken * if fan is well, to check how long time it became overheating after start, if it will be hotter and hotter * Call service provider and inform the details
Continuous alarm and concentrator not operating with power switch 'on'	* Unit overheating due to blocked air intake *Safety valve of compressor start self-protect system	Turn off the power switch, to re-open till the ball of flow meter drops . Restart the concentrator If the issue repeats, call service provider
Power switch on, the equipment works well, but no oxygen out and ball of flow meter at the bottom not rising up	* the flow knob is too tight or broken * the equipment is leaking badly	* power switch on, turn the flow knob anticlockwise, to check the adapter whether move with the knob together, if not, means the flow meter is broken. Otherwise, check the leaking of this product * Call service provider and inform the details

Nebulizing too slow	* not assembled well or broken * outlet is not tighten * flow meter not turn to "0"	* re-assemble or replace nebulizer * surely tight the outlet * * adjust the flow meter to "0"
All other problems		Call service provider

## 9. QUALITY WARRANTY

MODEL	WARRANTY PERIOD
OC-E80	3 Years
OC-E100	3 Years

EN

All after-sale service commitment shall be fulfilled by the supplier, dealer or the appointed maintenance company.

During the warranty period, if the product fails under conditions of normal use, the damages to the equipments are not caused by man-made reasons, the supplier, dealer or the appointed maintenance company is responsible to the maintenance or replacement free of charge.

The warranty does not include consumables like oxygen tubing, filters etc. The warranty does not cover the breakdown or damages caused by improper operation, abuse of the product, accident, or to products damaged by reason of repairs made to any components without the specific consent of the supplier. The users should not dismantle the equipment by themselves in the warranty period.

This warranty does not include normal wear and tear or shipping charges. SysMed and the dealer's sole obligation and your exclusive remedy under this warranty shall be limited to such repair or replacement. SysMed shall not be liable for any consequential or incidental damages whatsoever.


After the guarantee period, or damages excluded in the warranty, the supplier, dealer or the appointed maintenance company will supply customers the paid service.


# ATTACHMENT I CIRCUIT DIAGRAM

## OC-S (220V)

EN

## OC-E, OC-P (220V)

## 10. GLOSSARY

	Consult instructions for use		Follow Instruction for Use
	Type BF Applied Part		Temperature limit
	Power on		Power off
	Caution		Fragile, handle with care
	Stacking Limit By Number		Manufacturer
 No Open Flames	No open flames		Keep dry
 No Smoking	No smoking		This Way Up
IP21	Drip Proof Equipment		AC Power
 0123	CE certification mark		Suitable for homecare use

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